# pco.dicam C1 intensified 16 bit scMos camera









### pco.dicam C1

After more than 30 years of experience with image intensified cameras, we are proud to introduce the new pco.dicam C1 to you. The pco. dicam C1 is the first intensified camera system which exploits the full performance inherent to **scientific CMOS** sensor technology.

It is the optical coupling of 25 mm high resolution image intensifiers with an outstanding high efficiency tandem lens system to a 16 bit 4 Mpixel **sCMOS** sensor which makes the camera so unique. The 10G fibre optic based data interface guarantees you uncompressed and robust 16 bit data transfer of 106 full frames per second via optical fiber over virtually any distance.




106 frames/s @ full 4.2 Mpixel resolution	high frame rates at high resolution for imaging of dynamic processes
> 7000 frames/s @	kHz scan rates for spectroscopic applications
reduced resolution	KITZ SCALLTATES for spectroscopic applications
1.1 e- readout noise	lowest readout noise of any gated intensified camera system
1.1 e- readout noise	lowest readout hoise of any gated intensilled carriera system
16 bit digitization	taking advantage of the higher dynamic range possible from high end image intensifiers
25 mm high resolution image intensifier	doubles the optical resolution of conventional 18 mm image intensifiers
optical coupling via ultra-speed tandem lens	outstanding image quality with high transmission efficiency and no artifacts
tandem lens with 0.53 : 1 image scaling	full 25 mm diameter of intensifier output is lossless imaged onto sCMOS sensor
10G fibre optic based data interface	fiber optical interface virtually covers any distance without deploying additional interface converters or signal amplifiers with immunity to EMI
880 MByte/s image data rate	highest sustained image data rate of any intensified camera system on the market; no limitations for recording duration
double image mode with	two consecutive full resolution images with a configurable
300 ns interframing time	minimum interframing time of 300 ns
4.2 MPix sCMOS sensor	overcomes CCD limitations in terms of speed and sensitivity
enhanced	fast MCP gating for improved extinction ratio for the blue
extinction ratio gating	and uv part of the spectrum
additional optical trigger input	robust trigger transmission over long distance in EMC critical environments
EF lens control -optional-	convenient remote lens control for camera systems inaccessible during an experiment
selected highly homogeneous image intensifiers	integrated best image intensifier quality available on the market
< 50 ns trigger to exposure start delay	ultra fast camera reaction to trigger event
4 ns gating with 25 mm intensifier	captures fast transient phenomena
extensive and highly precise IN/OUT signaling	allows for perfect synchronization in any experimental set-up as timing master or slave
configurable delay in steps of 1 ns	flexible adaptation to synchronization needs



### pco.dicam C1

### >> applications

laser induced incandescence (LII) | shock wave physics | laser induced breakdown spectroscopy (LIBS) | particle image velocimetry (PIV) | time resolved spectroscopy | plasmaphysics | laser induced fluorescence (LIF) | ballistics | combustion




10G fibre optic based interface 8

# pco.dicam C1



### » sCMOS image sensor


type of sensor	SCIENTIFIC CIVIOS (SCIVIOS)	
resolution (h x v)	2048 x 2048 active pixel	
pixel size (h x v)	6.5 µm x 6.5 µm	
sensor format / diagonal	13.3 mm x 13.3 mm / 18.8 mm	
shutter mode	single image	
	double image	
MTF <sup>1</sup>	76.9 lp/mm (theoretical)	
fullwell capacity	15 000 e- for P46 phosphor	
	30 000 e- for P43 phosphor	
readout noise <sup>2</sup>	1.1 med / 1.5 mms e- single image	
	2.2 med / 2.5 ms e- double image	
dynamic range	13 600 : 1 (82.7 dB) for P46 phosphor	
	27 200 : 1 (88.7 dB) for P43 phosphor	
quantum efficiency	58 % for P43 peak emission @ 545 nm	
	57 % for P46 peak emission @ 530 nm	
spectral range	300 nm 1000 nm	
dark current <sup>3</sup>	< 0.6 e <sup>-</sup> /pixel/s @ 7 °C	
DSNU	1.0 ms e <sup>-</sup>	
PRNU	< 0.6 %	
anti blooming factor	1:10 000	



### >> frame rate table<sup>4</sup>


### ) quantum efficiency


2048 x 2048	106 fps
2048 x 1024	210 fps
2048 x 512	414 fps
2048 x 256	807 fps
2048 x 128	1535 fps
2048 x 64	2795 fps
2048 x 32	4739 fps
2048 x 16	7266 fps
1920 x 1080	199 fps
1600 x 1200	180 fps
1280 x 1024	210 fps
640 x 480	441 fps
320 x 240	858 fps



 $<sup>{1\</sup>over 2}\,{\rm Modulation\,transfer\,function}$   ${2\over 2}\,{\rm The\,readout\,noise\,values\,are\,given\,as\,median\,(med)\,and\,root\,mean\,square\,(rms)\,values\,due}$ to the different noise models, which can be used for evaluation. All values are raw data without  $\frac{1}{2}$ 

any filtering.  $^3$  Measurements with dark current compensation  $^4$  Exposure time < 1  $\mu s$ 






# pco.dicam C1

### >> camera system

frame rate	up to 106 fps
	@ 2048 x 2048 pixel
dynamic range A/D <sup>5</sup>	16 bit
pixel scan rate	286.0 MHz
binning horizontal	x1, x2, x4
binning vertical	x1, x2, x4
region of interest (ROI)	horizontal: steps of 4 pixels vertical: steps of 1 pixel
non-linearity	< 1 %
cooling method	+7 °C stabilized, 1 stage peltier with forced air (fan)
input signals	optical trigger (FOL), electrical trigger, arm input (TTL level, BNC connectors), gate disable (high-speed TTL input, BNC connectors)
output signals	gate/expos out monitor, user monitor output (TTL level, BNC connectors)
time stamp	in image (1 µs resolution)



### single image mode

exposure times	4, 10 ns fixed, 20 ns 250 ns (1 ns steps), 250 ns 1 s (10 ns steps)
delay times	0 ns 250 ns (1 ns steps), 250 ns 1 s (10 ns steps)
maximum repetitionwith external gating	200 kHz sustained, 3.3 MHz burst
insertion delay	
trigger input to exposure out	19 ns
trigger input to optical open	49 ns
jitter	
trigger input to exposure out	35 ps rms
trigger input to optical open	150 ps rms

### double image mode

exposure times	20 ns 1 ms (in 10 ns steps)
delay settings	0 ns 10 ms (in 10 ns steps)
interframing time	300 ns 10 ms (in 10 ns steps)

### >> general

power supply	18 28 VDC
power consumption	35 40 W
weight	7 kg
operating temperature	+ 10 °C + 40 °C
operating humidity range	10 % 80 % (non-condensing)
storage temperature range	- 10 °C + 60 °C
optical interface	F-mount
	optional: C-mount, Canon mount
lens remote controller	electronic control for Canon EF lenses
CE / FCC certified	yes




<sup>&</sup>lt;sup>5</sup> The high dynamic signal is simultaneously converted at high and low gain by two 11 bit A/D converters and the two 11 bit values are sophistically merged into one 16 bit value.




# pco.dicam C1

### ) image intensifier

type	HighRes MCP <sup>6</sup> (6 µm channel)
input window	synthetic silica, borosillicate
	(MgF₂ optional)
photocathode material	S20, GaAs, GaAsP (others on request)
image intensifier	6 μm
pitch distance	
image intensifier MCP type	single stage low resistance MCP for
	high strip current
MCP operational modes	continuous
	gated for enhanced extinction ratio
image intensifier diameter	25 mm (18 mm optional on request)
phosphor screen material	P43, P46
output window	glass
image intensifier	> 50 lp/mm @ 5 % MTF typical
system resolution	(depends on phosphor)
shortest gating time	4 ns



### >> image intensifier input window




Typical transmittance of image intensifier input window materials


To make use of the good UV-sensitivity of S20 photocathode material the standard input window is made of synthetic silica for transmission down to 180 nm.

GaAs and GaAsP photocathodes are deposited on borosilicate glass.

### pco.dicam C1



### >> image intensifier photocathode characteristics



Spectral sensitivities of different photocathode materials: S20 (multialkali), GaAs, GaAsP

data courtesy of Hamamatsu Photonics

photocathode material	peak wavelength [nm]	typical quantum efficiency at peak wavelength [%]	dark counts [s <sup>-1</sup> /cm <sup>2</sup> ]
S20 (multialkali)	250	20	1500
GaAs	650	30	30 000
GaAsP	500	55	10 000

data courtesy of Hamamatsu Photonics

### image intensifier phosphor

phosphor	phosphor decay (typ.) to		peak	typical
priosprior	10 %	1 %	emission	efficiency
P43	1 ms	4 ms	545 nm	100 %
P46	0.2 - 0.4 µs	2 µs	530 nm	30 %

You can combine all photocathode materials with P43 or P46 phosphor. Whereas the P43 phosphor has a much brighter emission than the P46 phosphor, it has a rather long decay time, i.e. the time required till the phosphor emission fades out after the excitation by electron bombardement has been stopped. This decay time is therefore critical for fast image repetition rates primarily in double image application or when operating the camera in spectroscopic mode with line rates in the kHz range.


### pco.dicam C1



### » optical coupling lens system

"ultra-speed tandem lens" between image intensifier & sCMOS

transmission efficiency	> 30 %
vignetting	< 3 %
resolution	> 60 lp/mm
scaling rates	B=0.53 for 25 mm intensifier



The projected image circle is completely covered by 2048 x 2048 6.5 µm pixels of the sCMOS detector – cf. image left. There is no "waste" of valuable intensifier area. As a consequence the four corners of the sCMOS sensor remain black. For a fast scan of just a few vertically centered lines – the camera module allows you to achieve more than 7.000 fps for such a ROI - the full line length of 2048 pixels is available.





### >> camera interface

data transfer	Camera Link HS, FOL cable, frame grabber (Single F2,1X1, S10)
input signals	optical trigger (FOL), electrical trigger, arm input (TTL level, BNC connectors), gate disable (high-speed TTL input, BNC connectors)
output signals	gate/expos out monitor, user monitor output (TTL level, BNC connectors)

### software

With pco.camware you control camera settings as well as image acquisition. Also, you can chose various formats for your files (Windows 7 and later).

Visit www.pco.de for a free software development kit (SDK). It includes a dynamic link library for user customization and integration on PC platforms. Drivers for popular third party software packages are also available.

# pco.dicam C1


### >> lens remote controller

The optional Canon lens control adapter enables you to connect electronic EF- and EF-S Canon lenses allowing to remote control focus and aperture of those lenses.



### dimensions

F-mount and C-mount lens changeable adapter.



All dimensions are given in millimeter.

### >> camera view





#### customization

### pco.dicam C1

#### » possible combinations

photocathode	input window	phosphor	order code
S20 selected*	synthetic silica	P46	70108001003
		P43	70108001002
GaAs standard*	borosilicate	P46	70108001006
		P43	70108001004
GaAs selected*	borosilicate	P46	70108001007
		P43	70108001005
GaAsP standard*	borosilicate	P46	70108001010
		P43	70108001008
GaAsP selected*	borosilicate	P46	70108001011
		P43	70108001009

 $<sup>^{\</sup>star}$  Image intensifiers with GaAs and GaAsP photocathode are available in two quality grades:

Standard:

quality specified for central 16 mm x 16 mm square region corresponding to 1300 x 1300 pixel sCMOS sensor resolution

Selected:

quality specified for 24.9 mm diameter area corresponding to full 2048 x 2048 pixel sCMOS sensor resolution, extinction ratio 10 times higher than standard grade, image intensifiers with S20 photocathode exclusively come in selected grade quality, contact our technical sales

team for further details on the two quality grades

### >> select optical interface

C-mount
F-mount
EF lens control

### >> select FOL

type of data interface FOL module in camera and frame grabber

SM SFP+ up to 10 km

MM SFP+ up to 300 m

FOL cable length default: 10 m

### Need more help?

Get in touch with our experts, tell us what you want to achieve and let us help you find the best customization for your application!

### find us

### europe PCO AG

Donaupark 11 93309 Kelheim, Germany

+49 9441 2005 50 info@pco.de pco.de

#### america

PCO-TECH Inc. 6930 Metroplex Drive Romulus, Michigan 48174, USA

+1 248 276 8820 info@pco-tech.com pco-tech.com

#### asia

PCO Imaging Asia Pte. 3 Temasek Ave Centennial Tower, Level 34 Singapore, 039190

+65 6549 7054 info@pco-imaging.com pco-imaging.com

### china

Suzhou PCO Imaging Technology Co., Ltd. Room A10, 4th Floor, Building 4 Ascendas Xinsu Square, No. 5 Xinghan Street Suzhou Industrial Park, China 215021

+86 512 67634643 info@pco.cn pco.cn



















for application stories please visit our website



