



Twin*Guard*

-86°C Ultra-Low Temperature Freezers





575 L / 715 L

PHCbi's **Twin***Guard* -86°C chest freezers deliver the unsurpassed level of safety and sample security that is critical for valuable biological sample storage alongside exceptional ease of use and data monitoring.

Ultimate Sample Protection

The Dual Cooling System "TwinGuard" comprised of independent dual cooling circuits delivers -86°C, and one side of the cooling circuits maintains a chamber temperature at -65°C. The freezer minimises the risk of compromising valuable samples.

Improved Temperature Performance

TwinGuard chest freezers maintain better temperature stability and uniformity during routine door openings than equivalent upright models as less cold air escapes the freezer while the door is open.


Enhanced Use & Intelligent Security

A colour LCD touch panel allows full user control, even with gloved hands, while the USB port makes transferring logged data to a PC convenient. The filterless design of the freezers reduces routine maintenance.



Valuable Sample Storage

With two independent systems, each sustaining ultra-low temperatures, your samples remain safe. Optional CO_2 backup provides an extra layer of protection.



Predictive Maintenance

Performance monitoring ensures the freezer is operating within specifications. This protects the compressors and condenser fans. External stress on the compressors or fans, due to ambient temperature abnormality or excess load, is displayed for immediate notification.



Restricted Access Laboratories

Significantly extended critical time to react to any sudden mechanical failure. Ideal for high security applications.

Twin Guard

C Ultra-Low Temperature Freezers



Filterless Design


The filterless construction of the freezers reduces routine maintenance time by eliminating the need for regular cleaning of filters.

Improved Temperature Performance


The TwinGuard -86°C chest freezer design provides optimal temperature stability and uniformity during routine lid openings, minimizing the loss of cold air even when the freezer lid is open.

Performance Data


AT30°C Pull-down & Pull-up Temperature



Dimensions



MDF-DC500VX



MDF-DC700VX

Temperature control range 2) °C -50 to -86 Control Controller Microprocessor, non-volatile memory Display LCD Touch Screen Pt-1000 Temperature sensor Refrigeration Refrigeration system Independent Dual-Cooling 2 x 1100 Compressors Refrigerant HFC mixed PUF / VIP PLUS Construction Exterior material Painted Steel Interior material Stainless Steel Outer doors/lid qty Outer door lock Inner doors/lid 3 pieces (Styrofoam) Access port Access port position Access port diameter Ømm 17 (V = Visual Alarm, B = Buzzer Alarm, R = Remote Alarm) V-B-R High temperature V-B-R Low temperature Filterless design Door open V-B MDF-DC500VX-PE MDF-DC700VX-PE MDF-DC700VXC-PA **Electrical and Noise Level** Power supply Frequency Hz 50 Noise level dB (A) 52 Options Liquid CO₂ back-up MDF-UB5-PW MTR-G85C-PE — Chart paper: RP-G85-PW — Ink pen: PG-R-PW - Circular type MTR-85H-PW — Chart paper: RP-85-PW — Ink pen: DF-38FP-PW - Continuous strip type - Recorder housing: MDF-S3085-PW **Optional Communication Systems** Ethernet interface (LAN) 41 MTR-I 03-PW Digital interface (RS232C/RS485) 4 MTR-480-PW **Quality Management System**

- 1) External dimensions of main cabinet only, excluding handle and other external projections - See dimensions drawings on website for full details.
- ^{2]} Air temperature measured at freezer centre, ambient temperature +30°C, no load.
- ^{3]} Nominal value Background noise 20dB[A]

Certification

External dimensions (W x D x H) 1

Internal dimensions (W x D x H)

Cooling performance 2

Temperature setting range

Net weight

Capacity

4) Only for MTR-5000 (data acquisition system) users.

Caution: PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for any loss or damage to the contents stored in the

IS09001

MDF-DC500VX-PE | MDF-DC700VX-PE | MDF-DC700VXC-PA

358

-86

-50 to -90

2010 x 845 x 1070

1190 x 640 x 756

575

328

416

litres

kg

" boxes

°C

2300 x 845 x 1070

1480 x 640 x 756

520

366



Preservation Equipment, Experimental Environment Equipment, Dispensary Equipment, Culturing Equipment and Drying & Sterilising Equipment for General Laboratory use

Unit: mm

The management of the design, development, production and servicing of the above. PHC Corporation, Biomedical Division

1-1-1 Sakada, Oizumi-machi, Ora-gun, Gunma 370-0596, Japan





PHC Corporation Biomedical Division is certified for:

60

Environmental management system: IS014001

DISTRIBUTED BY:



電話:02-2980-7936 傳真:02-8988-3071 地址:新北市三重區五華街1巷26號1樓

信箱:service@sagevision.com.tw 回题回

網站:http://sagevision.com.tw/ 高粱



PHC Corporation